

LM-MultiSeptor

Effective and easy-to-use saliva ejector

LM-MultiSeptor keeps the area being treated dry – both lingual and buccal. It functions equally well on both the upper and lower jaw.

- can easily be formed to the desired shape
- comfortable for the patient
- leaves ample space for procedures
- an invaluable aid when working unassisted.

The MultiSeptor is compatible with all standard suction units.

Care: disinfect, rinse, remove connections clean, and sterilize. Connections must be cold-sterilized.

LM-MultiSeptor, 40 cm

LM 2020

LM-MultiSeptor, 50 cm

LM 2020A

Spare Parts

Replacement tubes = 5 tubes (40 cm) + spring
LM 2021

Replacement tubes = 5 tubes (50 cm) + spring
LM 2021A



Brush Application Tips

Brush Application Tips are for one time use only!

- used for applying liquid treatment agents to all dental surfaces or cavities
- excellent for placing cavity liners, bonding and priming agents and etchants
- by bending the shank of the application tip the access into the cavities can be optimised

Application tips are available in three different colours (red, blue, yellow) to ease the identification and to coordinate multi-step procedures.



Brush Application Tip, 50 pieces

LM 758Y yellow

LM 758B blue

LM 758R red

Use Brush Application Tips together with LM-BrushHolder (LM 481-758). See page 38 for more information on the BrushHolder.

Code Rings

- LM 5000 blue
- LM 5001 green
- LM 5002 orange
- LM 5003 red
- LM 5004 yellow
- LM 5005 black
- LM 5006 white
- LM 5007 brown
- LM 5008 lavender
- package 50 pcs/colour



Convenient LM-Servo cassettes

Save time, money and personnel

- made of an exceptionally durable special plastic
- cassettes tolerate all common detergents, washers, and methods
- max. 180°C / 355°F
- standard dimensions – can be kept in all commonly used instrument cabinets, instrument maintenance equipment, and covered trays
- one-piece, light-weight

LM-Servo 5 LM 6650

- for five instruments
- dimensions: 180 x 86 x 28 mm



PAT.

LM-Servo 8 LM 6680

- for eight instruments
- dimensions: 180 x 136 x 28 mm



PAT.

Code buttons

On both of the cassette's long sides there are holes for three code buttons

- buttons sold separately
- package 20 pcs/colour



LM-ProColor instrument disinfectant

Effective protection

LM-ProColor is a new aldehyde-free disinfectant for hand instruments that combines a wide effective spectrum, high cleaning power and composition that protects various materials. LM-ProColor is also suitable for use in ultrasonic cleaners.

Safe and long-lasting

LM-ProColor is a long-lasting concentrate, whose pleasant aroma and colour enhance user comfort. One two-liter bottle produces 100 liters of cleaning solution.

Composition: 100 g LM-ProColor contains: 9 g alkyl-benzyl-dimethylammonium chloride, 0.1 g guanidine compound, non-ionic surfactants, alkaline detergents, corrosion inhibitors, adjuvants in aqueous solution.

Effective spectrum: bactericidal inclusive of TB bacteria, fungicidal, virucidal (enveloped viruses, HBV¹⁾, HCV²⁾, HIV³⁾).

1) by electron microscopy
2) Test for BVDV surrogate virus
3) Effectiveness against HBV implies effectiveness against HIV.

LM-ProColor LM 4100

2 l concentrate

