

The image shows several orthodontic instruments with colorful handles (orange, blue, yellow, blue, green) and metal shafts. One instrument has a curved mirror handle, another has a hook-like tip, and a third has a small blue tip. A pair of glasses is visible in the top right corner. A yellow banner at the top contains the title 'Orthodontie'.

Orthodontie

- Ortho-Spiegelhandgriff
- Ligaturenadapter
- Drahtbiegeinstrument

Orthodontie

**Ortho-Spiegelhandgriff
LM 25-26 XSi**

- nur mit dem Handgriff LM-ErgoMax XSi erhältlich
- zum Messen von horizontalem und vertikalem Überbiss
- 25 mm lange Millimeterskala




**Ligaturenadapter
LM 414-415**



**Ligaturenadapter - Scaler U15
LM 414-150**




**Gummiligaturenadapter
LM 416-417**



**Drahtbiegeinstrument
LM 412-413 XSi**

- nur mit dem Handgriff LM-ErgoMax XSi erhältlich



Das Ende eines Drahtbiegeinstruments abzubiegen kann recht problematisch sein. Die Platzverhältnisse sind beengt und der Draht ist starr. Das doppelendige LM-Drahtbiegeinstrument erleichtert das Biegen. Die Enden in Form abgeflachter Kugeln verfügen über 1 mm große Öffnungen, durch die sich der Draht problemlos einführen und dann biegen lässt. Die Klinge ist lang und erreicht problemlos selbst schwerzugängliche Stellen.