

The image features a collection of orthodontic instruments. At the top, there is a red handle, a blue handle, and a pair of glasses. Below these, a yellow handle is visible. In the lower section, there are blue and green handles, and a pair of pliers. The tools are arranged diagonally across the frame. A yellow banner at the top contains the title 'Orthodontie'. A list of three items is positioned in the middle-right area.


Orthodontie

- Manche de miroir, ortho
- Presse-ligature
- Presse-arc


Orthodontie

**Manche de miroir, ortho
LM 25-26 XSi**

- disponible uniquement avec le manche LM-ErgoMax (XSi)
- détermination de l'occlusion excessive dans le sens horizontal et vertical
- graduation millimétrique de 25 mm de long

**Presse-ligature
LM 414-415****Presse-ligature - Scaler U15
LM 414-150****Applicateur de ligature en caoutchouc
LM 416-417****Arch Wire Tucker
LM 412-413 XSi**

- disponible uniquement avec le manche LM-ErgoMax (XSi)



L'inclinaison de l'extrémité d'un arc peut être problématique. L'espace de travail est limité et l'arc rigide. LM-Arch Wire Tucker facilite l'inclinaison par sa double extrémité. Sa tête ronde aplatie est dotée d'orifices de 1 millimètre dans lesquels le fil peut être aisément introduit et tourné. La lame est allongée et s'adapte parfaitement aux espaces les plus serrés.